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IPG Wireless Device Alert Conditions Hide Alert Group Conditions

Note: Any alert condition group setting may be overridden for an individual patient. Customized settings are assigned within the Patient Profile.

Display descriptions of alert conditions

	Red Alert	Yellow Alert	Website-only Alert
Clinical Management Alerts			
AT/AF Daily Burden > Threshold	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Average Ventricular Rate during AT/AF	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Monitored VT Episode Detected	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Cumulative Right Ventricular Pacing > 40%	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Total Ventricular Pacing < 90%	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
OptiVol Fluid Alert	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Lead/Device Integrity Alerts			
Low Battery Voltage Recommended Replacement Time	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Atrial Pacing Impedance Out of Range	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Right Ventricular Pacing Impedance Out of Range	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Left Ventricular Pacing Impedance Out of Range	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
2nd Left Ventricular Pacing Impedance Out of Range	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Atrial Capture Management High Threshold	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Right Ventricular Capture Management High Threshold	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Left Ventricular Capture Management High Threshold	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
2nd Left Ventricular Capture Management High Threshold	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Non-Programmable Alerts			
Device Reset *	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Diagnostic Reset *	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Note: Alerts identified with an asterisk () are non-programmable alerts.

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Medtronic CareAlert Setup

Lead/Device Integrity Alerts	Device Tone Enable - Urgency	Patient Home Monitor
RV Lead...	2 of 2 On - High	Yes
Lead Impedance Out of Range...	4 of 4 On - High	2 of 2 On
Low Battery Voltage RRT...	On-High	4 of 4 On
Excessive Charge Time EOS...	On-High	On
VF Detection OFF, 3+ VF or 3+ FVT Rx Off.	On-High	On

Alert Time... 08:00

Clinical Management Alerts...

Demonstrate Tones... **Undo Pending** **OK**

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
In-clinic follow-up

Patients make an appointment.

Remote follow-up

Patients are scheduled for their remote transmission through the clinic's scheduling system or CareLink.

Medtronic CARELINK® NETWORK Help | Resources | View Profile | Sign Out

HOME TRANSMISSIONS **MANAGE MY PATIENTS** MANAGE MY CLINIC Clinic: DEMOZEDU

[All Patients](#) [New Patients](#) [Add a Patient](#) [Pending Patient Transfers](#) [Transfer In](#) [Advanced Search](#) [Transmission Schedule](#) [Equipment Orders](#)

Manage My Patients

Schedule These Patients

To schedule the following patients for a common transmission date, select a date and click "Schedule". Uncheck the box for any patient you do not wish to schedule at this time.

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

Medtronic, Medtronic logo, and Further, Together are trademarks of Medtronic. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Medtronic is under license. ™Other third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

medtronic.eu

UC202006339 PO ©2020
Medtronic. All Rights Reserved.

Medtronic