



HeartLogic™ Heart Failure Diagnostic

In-Service

HeartLogic™
Heart Failure Diagnostic

Heart failure is a major clinical concern that impacts patients' quality of life and healthcare economics

Overview

HeartLogic™

MultiSENSE Study

Case Studies

Heart failure hospitalizations are associated with an **8-fold increase in mortality** and a **9-fold increase in recurrent hospitalizations**¹




Figure adapted from Goodlin¹

HF hospitalizations are the **biggest outflow** of money

- Accounting \approx 60% of healthcare expenditure²

Heart failure symptoms of decompensation are possible to identify

Overview

HeartLogic™

MultiSENSE Study

Case Studies

Symptoms can vary in presentation and time course, nevertheless, several patterns may precede a decompensation

“Typical” Progression of Heart Failure Decompensation

Cardiac output ↓

Left atrial pressure ↑

Pulmonary Fluid ↑

Symptoms ↑
Activity ↓

Sympathetic nervous activity ↑

Weight ↑

Hospitalization

Heart Sounds

HRV

Thoracic impedance
Respiration
Sleep incline

Weight

Activity

Boston Scientific sensors enable the early detection of heart failure decompensation symptoms



Overview

HeartLogic™








MultiSENSE Study

Case Studies

Sensors are intended to represent typical in-office tests and questions

Our sensors:

What Clinicians ask/do during a physical exam

 Heart Sounds _____	Listen to the heart S3 heart sound
 Heart Sounds _____	Listen to the heart S1 heart sound
 Thoracic Impedance _____	Take chest X-ray for signs of pulmonary edema
 Respiratory Rate _____	"Are you out of breath? Have difficulty breathing?"
 Activity Level _____	"Are you able to get your mail/go upstairs?"
 Weight _____	"Have you gained weight?" (check leg or abdominal swelling)
 Night Heart Rate _____	Is resting heart rate elevated?

Introducing HeartLogic™ Heart Failure Diagnostic

Boston
Scientific

Overview


HeartLogic™

MultiSENSE Study

Case Studies

HeartLogic™ enables **proactive care** of heart failure patient management **with a multifactorial approach**

Incorporates Multiple Sensors
with a Single Composite Alert



Available on LATITUDE NXT for
patients with **Resonate™** family of ICDs
& CRT-Ds



The MultiSENSE study results¹ demonstrated:
High sensitivity of 70 % for detecting heart failure events
Weeks of advance notice of a potential heart failure event
Low burden of less than 2 alerts per patient per year.

HeartLogic™
Heart Failure Diagnostic

1) Boehmer, J et al., JACC-HF, 2017;5(3),2 1 6 – 2 5.
CRM-495201-AB

Benefit of Multifactorial Approach

Overview



HeartLogic™

MultiSENSE Study

Case Studies

Patient A — Two Observed Cases — Patient B

Which patient
had a Heart Failure Event?



Benefit of Multifactorial Approach

Overview

HeartLogic™



MultiSENSE Study

Case Studies

Patient A — Two Observed Cases — Patient B

Multi-sensor Changes before a **HF Event**

Impedance-only Change with **NO** Event



HeartLogic Uses Multiple Sensors

HeartLogic™ Heart Failure Diagnostic

Overview

HeartLogic™

MultiSENSE Study

Case Studies

HeartLogic was validated¹ to detect the early warning signs of worsening heart failure by combining data from **5 sensors** into a single composite index.

Multiple Sensor Measurements



Heart Sounds
S1 & S3



Impedance
Thoracic



Respiration
Rate & Volume





Activity
Time Spent Active



Heart Rate
Night

Combined into a single, simple index with alert

HeartLogic™ Heart Failure Index



HeartLogic™ Includes:

- a** Composite HeartLogic™ Index trend
- b** Actionable HeartLogic™ Alert
- c** Configurable HeartLogic™ Threshold
- d** Heart Failure Management Report with HeartLogic™ Data

Heart Sounds; new to Resonate Platform, Unique to Boston Scientific

Overview


HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

Heart Sounds were traditionally evaluated based on physician audibility using a stethoscope



- In a normal subject, the first heart sound S1 happens shortly after the Q-wave in an EKG, produced by the closing of atrioventricular valves.
- S3 is caused by rapid filling of blood against a stiffer ventricular wall during passing filling. It may not be audible in normal subject

Heart Sounds; new to Resonate Platform, Unique to Boston Scientific

Overview



HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

Heart sounds are measured from the accelerometer and reveals signs of elevated filling pressure and weakened ventricular contraction via S3 and S1 heart sounds.



Worsening heart failure may be associated with an increase in S3 or a decrease in S1, or both

Device-measured S3 more efficient than auscultation

Overview



HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

Device-based **S3** includes audible and sub-audible frequencies and therefore can **pick up vibrations that cannot be detected by the human ear**



- Device-based S3 was significantly louder in patients with HF events

- Auscultated S3 was mostly absent and no different between groups

Device-based S1 and S3 measured during clinically stable periods

Overview



HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

S3 and S1 amplitudes were **significantly different** in **pre-event windows** from those during clinically stable periods in ambulatory heart failure patients



A clinically stable period is defined as a period between office visits when: a) NYHA classification was unchanged, b) Weight change \leq 5 lbs (2.27kg), c) No adverse events were reported between the visits

- Presented at HFSA 2013; data from PRE-SENSE study, E. Klodas, et. al. Journal of Cardiac Failure, 2013, 19(8): S67
- Gardner RS, Nair DG, Eldadah Z, Stancak B, Thakur P, Wariar R, Zhang Y, An Q, Averina V, Sweeney R, Boehmer JP. Heart failure diagnostic sensor measurements during clinically stable epochs in ambulatory heart failure patients. Eur J of Heart Failure 2017;19(Suppl. S1):277(P1132). Cao M, Schulze C, Gardner R, An Q, Thakur P, Thompson J, Boehmer J. Device-measured third heart sound predicts heart failure events better than auscultated third heart sound. Europace 2017; 19 (suppl_3): iii332-iii333.
- Sriratanasathavorn C, Fisher J, Bransford P, Chan R, Sweeney R, Ahmed R, Averina V, An Q, Thakur P, Zhang Y, Boehmer J. Heart Failure Diagnostic Sensor Measurements Change Prior to Heart Failure Decompensation Events. Accepted at APHRS 2017
- TSUGUYA SAKAMOTO, REIZO KUSUKAWA, DONALD M. MACCANON, ALDO A. LUISADA and Ivan Harvey. Hemodynamic Determinants of the Amplitude of the First Heart Sound. Circulation Research. 1965;16:45-57, originally published January 1, 1965. <https://doi.org/10.1161/01.RES.16.1.45>. Thakur, P. H., An, Q., Swanson, L., Zhang, Y., and Gardner, R. S. (2017) Haemodynamic monitoring of cardiac status using heart sounds from an implanted cardiac device. ESC Heart Failure, doi: 10.1002/ehf2.12171.

Device-measured 3S HFE prediction Vs. auscultated third heart sound

Overview



HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

Device-measured **S3** has better correlation with the **risk of heart failure** events over one year than auscultated S3



This superiority may be attributable to the subjective nature of auscultation, limitations of the frequency range of human hearing, as well as limitations of a single snap-shot auscultation assessment.

Device-measured S3 significantly correlated to echo parameters

Overview

HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

In PRE-SENSE, a related feasibility study conducted on CRT-D patients:

- S3 amplitude was significantly correlated with EDT, EDR, E/A ratio and E/E' ratio, which are often used to assess LV diastolic function
- There was significant separation across almost all echo parameters between the low and high S3 groups (high = top quartile).

Measures	Low S3 (n=48)	High S3 (n=14)	p-value	
LVEF (%)	41.5	28.4	0.004	
LV Volume	LVEDD (cm)	5.57	6.71	<0.001
	LVEDV (mL)	134.8	182.6	0.005
	LVESD (cm)	4.35	5.81	<0.001
	LVESV (mL)	83.5	134.6	0.003
LA volume	LAEDV (mL)	72.1	101.1	0.006
Diastolic Filling (E and A waves)	E (cm/s)	62.8	89.0	0.006
	EDT (ms)	279.5	193.3	<0.001
	EDR (m/s ²)	2.63	5.02	<0.001
	E/A ratio	0.89	1.49	0.020
	E/E' septal ratio	13.2	19.4	0.003
E' septal (cm/s)	5.31	4.29	0.162	

HeartLogic Sensors

Overview


HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● **Respiration** ● Night Heart Rate ● Activity ● Thoracic Impedance ● Report

- The Respiratory Sensor uses transthoracic impedance measurements to collect respiration-related data, specifically respiration rate and tidal volume.
- Respiration rate and tidal volume are combined to produce additional respiration metrics.
- **Worsening heart failure may be associated with an elevated RR, an increase in RR, or an increase in day-to-day RR variability.**



Overview

HeartLogic™

MultiSENSE Study


Case Studies

● Heart Sounds ● Respiration ● **Night Heart Rate** ● Activity ● Thoracic Impedance ● Report

Worsening heart failure may be associated with **increases in Night Heart Rate**

Night Heart Rate

- Average Heart Rate from 12 am to 6 am
- Tends to capture close to resting heart rate for most patients
- Resting HR are known to be prognostic for cardiovascular outcomes
- Look for increases as a sign of heart failure worsening



Overview

HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● **Activity** ● Thoracic Impedance ● Report

Worsening heart failure may be associated with a **decrease in a patients activity level**

Activity

- Numbers of hours per day a patient is active.
- Overall status and fatigue.



Overview

HeartLogic™


MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● **Thoracic Impedance** ● Report

Thoracic Impedance²

- The Thoracic Impedance sensor measures impedance between electrodes on the RV lead and the pulse generator (PG) can.
- Thoracic Impedance may be associated with the fluid level in the patient's chest and may track fluid level changes.



Worsening heart failure may cause a patient's fluid level to increase, therefore, decrease in thoracic impedance. This information can be valuable when used in conjunction with a multifactorial approach to assess and treat HF patients¹

- Wang L, Lahtinen S et al. Feasibility of using an implantable system to measure thoracic congestion in chronic heart failure. Pacing Clin Electrophysiol. 2005;28:404-411
- Sriratanasathavorn C, Fisher J, Bransford P, Chan R, Sweeney R, Ahmed R, Averina V, An Q, Thakur P, Zhang Y, Boehmer J. Heart Failure Diagnostic Sensor Measurements Change Prior to Heart Failure Decompensation Events. Accepted at APHRS 2017.
- Yu CM et al. Intrathoracic impedance monitoring in patients with heart failure: correlation with fluid status and feasibility of early warning preceding hospitalization. Circulation, 2005 Aug; 9; 112(6): p841-8.
CRM-495201-AB

HeartLogic™ : Heart Failure Management Report

Boston
Scientific

Overview

HeartLogic™

MultiSENSE Study


Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● **Report**

HeartLogic collects information through the implanted device **without patient involvement** and provides data via **LATITUDE NXT 24/7**

HeartLogic™ Includes:

- a** Composite
HeartLogic™ Index trend
- b** Actionable
HeartLogic™ Alert
- c** Configurable
HeartLogic™ Threshold
- d** Heart Failure
Management
Report with
HeartLogic™ Data



High performing composite indicator for detecting worsening of heart failure using multiple physiologic measurements

HeartLogic™
Heart Failure Diagnostic

CRM-495201-AB
Boehmer, J et al., JACC-HF, 2017;5(3),2 1 6 – 2 5

Enabling HeartLogic™ :

Quick review

Overview

HeartLogic™

MultiSENSE Study

Case Studies

● Heart Sounds ● Respiration ● Night Heart Rate ● Activity ● Thoracic Impedance ● **Report**

There are essentially two steps required to enable HeartLogic:

- Program the device to ensure that Heart Failure Sensor Suite is activated. Heart Failure Sensor Suite is on by default
- Ensure that the patient with the implanted device is added to the Latitude™ NXT system, provided with a Communicator, and that the Latitude™ NXT system is programmed appropriately.

It will take up to 37 days from the time of implant to the point where it begins compiling Index numbers, allowing potential transmission of HeartLogic™ Alerts over Latitude™ NXT, assuming it is programmed ON.

In order to see HeartLogic data, remember that patients need to be on Latitude™ NXT and HeartLogic has to be activated for the patient.

Primary Results from the MultiSENSE Study

Overview

HeartLogic™

MultiSENSE Study

Case Studies

A **Multi-Sensor Algorithm Predicts Heart Failure** Events in Patients with Implanted Devices

- International, multi-center, non-randomized, clinical study designed to **develop** and **prospectively evaluate** a multi-sensor index and alert for the **early detection of worsening heart failure**.
- **New sensors created** with enhanced components and **novel data collection** and **processing techniques**.

Key inclusion criteria

- HeartLogicAge 18 or above
- Currently implanted with a COGNIS CRT-D system
- NYHA Class II, III or IV within the last 6 months

Key exclusion criteria

- Documented as pacemaker dependent
- History of appropriate Tachy therapy 1 week prior to enrollment
- Likely to undergo lead or PG revision
- Subjects that have received a heart or lung transplant
- Receiving mechanical circulatory transplant
- A life expectancy of less than 12 months

Primary Results from the MultiSENSE Study

Overview

HeartLogic™

MultiSENSE Study

Case Studies

Independent clinical events committee (CEC) Adjudication:

Heart failure Events (HFE)	<p>Primary cause of event was worsening heart failure and</p> <ul style="list-style-type: none">• Is admitted for HF and receives an augmented HF regimen with oral or intravenous medications, or• Receives unscheduled intravenous decongestive therapy that does not involve formal in-patient hospital admission, regardless of the setting
True Position Alerts	<ul style="list-style-type: none">• Onset before a usable HFEs• Recovery no earlier than 30 days before usable HFEs
HF Related Alerts	<p>Same onset and recovery window but broader set of HF events:</p> <ul style="list-style-type: none">• hospitalizations with a secondary cause of HF,• outpatient visits with a primary cause of HF and augmented oral medication changes,• HFEs that did not meet sensor data availability criteria or occurred within 45 days of device conversion
Unexplained Alerts	<ul style="list-style-type: none">• All other alerts

Primary Results from the MultiSENSE Study


Overview

HeartLogic™

MultiSENSE Study

Case Studies

Data was used to develop individual physiologic sensor trends and a **multi-sensor composite alert for worsening heart failure**



Primary Results from the MultiSENSE Study


Overview

HeartLogic™

MultiSENSE Study


Case Studies

Predefined endpoints



- Endpoint 1: **Sensitivity** for detecting usable heart failure events **>40 %**
- Endpoint 2: **Unexplained alert rate** (UAR) per patient year **<2.0**

Alert development



- **HeartLogic threshold** is configurable to user's preference for sensitivity and specificity
- By increasing sensitivity the UAR decreases

Primary Results from the MultiSENSE Study

Overview


HeartLogic™

MultiSENSE Study

Case Studies

The MultiSENSE Study data validated that HeartLogic™:

- Had **high sensitivity** of 70 % in detecting heart failure events
- Had a very **low burden** of less than 2 alerts per patient per year
- May have allow several **weeks of advanced notice** to clinicians of a potential event



HeartLogic™ Index Trends


Overview

HeartLogic™

MultiSENSE Study

Case Studies

MultiSENSE Study results showed that it could have allowed **weeks of advanced notice** to clinicians to a potential heart failure event.



HeartLogic™ Index
in patients **with** Heart
Failure Events

HeartLogic™ Index
in patients **without** Heart
Failure Events

Compared to 3 month baseline, HeartLogic Index was **statistically higher** beginning 29 days prior to event.


Event Rate Ratio Presented as a LBCT at ESC-HF 2017

Overview

HeartLogic™

MultiSENSE Study

Case Studies



- At nominal HeartLogic Threshold = 16
 - HF Event Rate was **10x higher** when **IN a HeartLogic Alert** versus **NOT in an alert**
- 17 % of patient-days IN Alert State

HeartLogic was better pronosticator than a baseline NT-proBNP


Overview

HeartLogic™


MultiSENSE Study

Case Studies

Event Rates as a Function of Daily HeartLogic Score



Event Rates as a Function of Baseline NT-proBNP




HeartLogic Significantly Augments The Prognostic Ability Of NT-proBNP Assessment

Overview

HeartLogic™

MultiSENSE Study

Case Studies



Half of patient follow-up in lowest risk group.

Event Rate Ratio compared to lowest risk group:

- LOW NT-proBNP, IN HL alert: **23.5**
- HIGH NT-proBNP, OUT HL alert: **8.0**
- HIGH NT-proBNP, IN HL alert: **50.0**

Conclusions

Overview

HeartLogic™

MultiSENSE Study

Case Studies

- HeartLogic was established and validated in the MultiSENSE trial with a **sensitivity of 70% and UAR of 1.47 alerts / patient-year at the nominal value of 16**
 - Other HeartLogic values can be tailored to the user's preference for sensitivity and specificity.
- **HeartLogic predicted risk** of HF event independent of baseline variables.
- Using HeartLogic threshold at 16, the median of detection was **34 days prior to a HFE**.
 - 89 % of true positive alerts occurred at least 2-weeks prior to a HFE.
- **Event Rate** was **10 times higher** when HeartLogic was **IN alert state** compared to OUT of alert state.
- Healthcare resources may be diverted to those with greater need.
- HeartLogic was a **better prognosticator than a baseline NT-proBNP**.

Thank You!

Boston
Scientific

HeartLogic™ shifts heart failure patient management from reactive treatment to **proactive care** that is validated in the MultiSENSE Study to have:

- **High sensitivity** of 70% for detecting heart failure events
- **Weeks of advance notice** of a potential heart failure event
- **Low burden** of less than 2 alerts per patient per year

HeartLogic™ incorporates multiple sensors with a single composite alert



Heart Sounds
S1 & S3



Impedance
Thoracic



Respiration
Rate & Volume



Activity
Time Spent Active



Heart Rate
Night

*multiple sensor
measurements
combined into
a single, simple
index with alert*

53

3 Mar 2017



Available on **LATITUDE™ NXT**
for patients with
Resonate™ family of
ICDs & CRT-Ds

